


Roll Boards for Lateral Transfer



Roll Boards work like a rolling conveyor belt to laterally transfer a person from bed to stretcher, or stretcher to imaging surface.

- Transfers are easily and safely done with 2 to 3 people, helping ensure no injuries to staff or patient
- Never leaves the area with the patient so it is always readily available for the next transfer
- Wipe down with any of your hospital-approved non-bleach disinfectant wipes
- Lightweight without any metal rollers
- Molded plastic handles for easy transport and for mounting the board on walls or back of doors
- 3-year manufacturer's warranty against any defects
- Disposable clear plastic covers are available for further protection in high fluid transmission environments



- Wall Mount Hook for Boards PN: 4086



- Standard or CT MRI Roll Board Disposable covers PN: 4010
- Surgi or All-Purpose Utility Board Disposable covers PN: 4014



- Dispenser and Hanger for Disposable Roll Board Covers (wall or door mount) PN: 4088



- Door Mount Hook for Boards PN: 4087

Also be sure to ask about our **SLIPP®** Lateral Transfer Solutions!

Slipp® is a unique patient moving device which makes transfers and repositioning easier than ever before on both the patient and caregiver! Unlike other patient moving devices, Slipp's unique design lets the patient (regardless of size) literally "slip" on a comfortable surface of sealed fluid that reduces normal friction to a minimum.

Better yet, there is no weight limit on any of the Slipp® units!!!

Easy, convenient, efficient and affordable, Slipp® is quickly becoming the transfer device of choice among leading hospitals, long term care facilities, and ambulance/EMS services nationwide.



DISTRIBUTED BY:

For more information, please contact:

Matt Center Phone: (603) 748-9762
matt@ergonomicsafetyinc.com

Roll Board types and suggested areas of use:



Standard Roll Board

- Size: 70"x20", fold in half to 35"x20" for easy transport and hanging of board
- Weight: only 6 lbs.
- Suggested Use Area: bed to stretcher transfers
PN: 4100



CT MRI Roll Board

- Size: 70"x16", fold in half to 35"x16" for easy transport and hanging of board
- Weight: only 6 lbs.
- Suggested Use Area: Stretcher to narrow imaging surface
PN: 4500



Surgi Roll Board

- Size: 48"x16", non-folding
- Weight: only 4 lbs.
- Suggested Use Area: OR table to stretcher
PN: 4400



Utility/All-Purpose Roll Board

- Size: 48"x20", non-folding
- Weight: only 4 lbs.
- Suggested Use Area: An economical alternative to having mixed sizes if a 3rd person is available to assist with the feet
PN: 4460

Our Roll Boards provide a safe and proven method to laterally transfer patients in any of your patient care areas.



No Weight Limit

DISTRIBUTED BY:

For more information, please contact:

Matt Center Phone: (603) 748-9762
matt@ergonomicsafetyinc.com